

# Smart Chill



Springs of Well-being


Skin Cooling System

# Smart Chill

Skin Cooling System

## Technical data


Description	Skin Cooling System
Treatment temperature	From 0°C up to 8°C (8°C is the typical used value)
Power Supply	230 Vac - 0.6 A - 50/60 Hz
Dimensions	18 cm (H), 20 cm (D), 31 cm (W)
Weigth	14.5 kg



[www.dekalaser.com](http://www.dekalaser.com)

### DEKA: THE SCIENCE OF WELL-BEING



Deka is the medical division of the Ei.En. Group, Italian leader in the production of laser equipment. Founded in 1981 it operates in more than 40 countries throughout the world. The Group's activity covers the majority of industrial applications and laser sources with powers of 6 kW. Deka produces laser equipment in compliance with specifications of Directive 93/42/EEC and adopts a strict quality assurance system certified by  according to the ISO 9001 and EN 46001 standards.



**DEKA M.E.L.A. srl**

Via Baldanzese 17 - 50041 Calenzano (FI) - Italy  
tel +39 055 8874942 - Fax +39 055 8832884  
[info@dekalaser.com](mailto:info@dekalaser.com)

*Product specifications are subject to change without notice.*  
COD:SC\_260204\_1